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Introduction

Compared with such instruments as are conventionally used in scoliosis surgery, internal fixators offer superior options for the segmental correction of positions, for fixation of the vertebral column through their repositioning levers which act on pedicle screws, and through their angle and rotation-secured principles of construction.

Based initially only on bi-segmental bracings, multi-segmental assemblies have long become a possibility.

For most fixators however, such longitudinal carriers as are required for definite immobilisation are already being used for repositioning. Most systems “have a reductive effect against the rod“. This would entail a few disadvantages in scoliosis surgery:

a) The correction is done in accordance with the pre-formed rod. Attainable targets are determined and consequently must be achieved.

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Introduction

b) The vertebral column is extended through the envisaged correction – de-rotation and elimination of lateral outward bending. In other words, during correction, the fastening points of the pedicle screws must be advanced on the longitudinal carrier, and during this process, a previously attained rotation correction can hardly be maintained.

c) The rod which is already lying in place during repositioning obstructs the step-wise release which is often required, as well as the adequate decortication of the laminae and joint obliteration for the preparation of spondylodesis.

As opposed to this, the krypton fixator has realised the principle of not using the in-laid rod as counter-bearing for the purpose of correction, but an instrument that is diverted outside the wound.


The process of the tri-dimensional correction and fixation of torsion scoliosis with the system krypton is illustrated as follows.

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Instruments

The correction instrument comprises the adaptors mounted on krypton pedicle screws, the repositioning spindle with a constrained terminal, and with one or several unconstrained spindle arms, the saddles for the absorption of the counter-bearing rods and the transversal tractions.

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In case of idiopathic scoliosis, the choice of the cranial or caudal end vertebral bodies of the instrumentation is made in accordance with the curved sections which are found as fixed in bending-absorptions (comparable with the popular surgical planning of VDS in accordance with Zielke).

Moreover, extended instrumentation may be required for other aetiologies, particularly in case of neurogenous scoliosis.

The selected end vertebral bodies are always instrumented on both sides, with krypton pedicle screws. Intermediate segments may well be skipped, above all, on the concave side. In the area of the end vertebral body, the pedicles on the convex side – the quantity of which is dependent upon the length of bend and the expected strength of the correction forces to be applied – are always instrumented.

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Operation technique

If no navigation system is applied, the x-ray image intensifier and the overturning of the operation table shall be aligned in such a way that the vertebra to be instrumented is presented precisely sideways in the horizontal radiation course of the C-arm. The pedicles on the concave side of the apex are often instrumented only after derotation because the deformation of the trunk, as caused by scoliosis, makes representation of the direction of the pedicle in advance difficult.

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The krypton pedicle screws are applied on both sides in the end vertebrae of the curve and on the convex side of the apex. Thereafter, the repositioning spindle is fastened with the spindle arms to the concave-side krypton pedicle screws of the end vertebrae. The end plate inclination of the end vertebrae can then be adjusted on the joints of the repositioning arms to meet the envisaged sagittal profile of the vertebral column after correction.

After occupying the convex-side pedicle-screws with transversal tractions, the actual correction can begin.

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The instrumented vertebra is turned segment-wise toward the concavity of the curve with the aid of the adaptor on the convex-side pedicle screw in the apex area. The position attained is secured with the open-end spanner, on each of the fixing joints of the respective transversal traction.

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The lateral outward curve of the bend which initially appears reinforced in this manoeuvre can be reduced through distance-reduction by the repositioning spindle. This eases the derotation process very much like the application of the T-wrench for adaptor for the extension of the repositioning arms on the krypton pedicle screws. Subsequently, the saddles are mounted on the joint-pieces of the repositioning arms as well as the counter-bearing rods. The hooks of the transversal tractions are attached to the counter-bearing rods.

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Operation technique

In the derotation of the apex area, the placement of the pedicle screws on the concave side is easier. These screws can be subsequently pulled on the lateral pass with the traction device for reduction scaffold of the krypton instrument. This reinforces the derotation and allows at the same time, for a controlled kyphosising of the sagittal profile in the area of instrumentation. Already at this point, the vertebral column is slightly extended and this extension is reduced again by the repositioning spindle.

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Correction step: Medialisation of the apex area

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At this point, the derotated apex area of the curve is centralised with the aid of the adjusting nuts of the transversal tractions which are attached to the counter-bearing rod. A concave-side release system makes this process easy and, at the same time, serves the purpose of preparing spondylodesis. The instrumented section of the vertebral column extends through the medial alignment of the apex of the curve, .

This is compensated for by a re-adjustment of the distance on the repositioning spindle.

Since the extension may lead to an expansion of the spinal cord, a wake-up test is required after this correction step as long as a continuous Spinal Cord Monitoring is not performed.

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Correction step: Setting the sagittal profile

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
To complete the correction process, the sagittal profile of the instrumented section of the vertebral column is finally set. In most cases this is realised through the respective manipulation of the krypton pedicle-screws attached to the repositioning spindle on the concave side, and subsequent renewed fixing of the joints on the repositioning arms (at the end vertebrae through the clearance on the saddle carrying the counter-bearing rods).

If necessary, an additional repositioning spindle can also be laid on the convex side. In rare cases of structural scoliosis, a partial facet resection is required on the convex side for the setting of an optimal lordosis in order to enable the adequate shortening of this range of joints.

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Permanent fixation and spondylodesis

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Prior to the assembly of the krypton rods, laminae and spinous processes are carefully decorticated and fine-crumbed autologous spongiosa is attached to the decartilaged facet articulation of vertebrae.

The krypton rods are prepared in accordance with the form and length found with the phantom rod. The rods are fitted with the necessary amount of krypton connectors and fastened to the krypton pedicle screws, whereby a centripetal advancement is more beneficial. The locking screws are positioned with the cardan screw driver. The rod is mounted on the convex-side which is occupied with more krypton pedicle screws.

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After the assembly of both longitudinal rods, further spongy bones are applied, one is centrally located and two terminal cross links are put in place. After very strong correction forces have been applied, it is advisable to have the two terminal cross links assembled prior to the removal of the repositioning instrument in order to reliably protect the convex-side krypton pedicle screws of the end vertebral bodies against lateral dislocation.

Early mobilisation and corset-free follow-up treatment is possible.