OP Technik  
*Surgical technique*

# Inhalt

## Einleitung

### Introduction

4

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### Indications and contraindications

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## Introduction

osmium is a cervical plate system for the anterior stabilization of the vertebral column. The implant immobilizes the fusion mass until the bone graft has consolidated.

The special feature of this system is the monocortical fixation with osmium screws. The osmium screw, consisting of a screw and a bolt, is expanded when the bolt is inserted.

The expansion mechanism makes superior stabilization and maximum anchorage possible, even if the bone is of reduced quality or if the same drill whole is used for a second screw placement in revision cases.

With the cancellous bone screws a bicortical fixation is achieved. The screws can be combined as per the surgeons preference.

osmium implants are made from a titanium alloy in accordance with the standard DIN EN ISO 9001.

The present surgical technique describes the working steps necessary for the application of the osmium. Also, the implants and the system instruments required for implantation are presented.

The surgical technique is not sufficient as the sole basis for the successful application of the osmium.

It is recommended to study and learn the operating technique with and from an experienced surgeon.

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The small pictures at the bottom of the page show the chronological application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-coloured.

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Indications

Contraindications

## Indications and contraindications

### Indications

osmium is applied in cases of cervical instabilities of different origin.


- Anterior disc resection
- Fractures
- Tumours
- Pseudarthrosis
- Revision surgery

### Contraindications

osmium must not be applied in cases of acute or chronic infection or in cases of proven intolerance to the materials used.


In case of poor medical or psychological general state of the patient which would be further aggravated through surgery careful consideration must be given.

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## ◀ osmium screw

The osmium screw (CS 1300-14, -16, -18) is a two-component system consisting of a screw and a bolt. The implant is available in different lengths (14, 16, 18 mm) with an outer diameter of 5 mm. The length is easily recognizable according to colour: 14 mm yellow, 16 mm green and 18 mm blue. Expansion of the screw is achieved by inserting the bolt. Hyperexpansion of the screw is not possible.

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## ◀ Cancellous bone screw

The cancellous bone screw (CS 1301-X) is used for the bicortical fixation. The outer diameter is 4 mm. The lengths of the screws are 12, 14, 16, 18, 20, 24 and 26 mm.

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## ◀ osmium plates

The patented plates are slotted at their cranial and caudal ends. The width of the plates is 20 mm.

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| Lochung | Länge                 |
|---------|-----------------------|
| 4-Loch  | 23, 27, 30, 33, 36 mm |
| 6-Loch  | 40, 44, 48, 52, 56 mm |
| 8-Loch  | 60, 66, 72 mm         |

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
▲ **Drill guides for repositor**  
(two are necessary)

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▲ **Tap for cancellous bone screws,** Ø 4.0 mm, for CS 1324

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
### ◀ Applying the repositor

The repositor (CS 1217-1) is positioned on the midline of the vertebral column and the drill guides (CS 1217-2) are inserted into the arms of the repositor.

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
### ◀ Drilling for reduction pins

In order to fix the reduction pins (CS 1210) a drill bit of Ø 2.0 mm (CS 1215) is used.

**Important:** Drill under C-arm control.

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### ◀ Placing of the reduction pins

After the removal of the drill guides (CS 1217-2) the reduction pins (CS 1210) are inserted into the arms of the repositor.

# Surgical technique – Monocortical instrumentation

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## ◀ Applying the reduction pins

The reduction pins are fixed by means of the screw driver (CS 1328).

**Important:** Introduce the reduction pins under C-arm control. After the reduction pins have been fixed, the intersomatic space may be distracted by means of the repositor (CS 1217-1 and CS 1217-3).

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## ◀ Bending of the plate

The osmium plate is contoured. Additional bending is possible if required by the patient's anatomy using the plate benders (CS 1327-1 and CS 1327-2).

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## ◀ Inserting of the plate


After distraction the intersomatic space can be cleared and a bone graft or cage can be inserted. In order to insert the osmium plate the repositor must be lifted slightly.

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### ◀ Compression

The achieved distraction may now be changed into compression with the lever arms for repositor (CS 1217-3).

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### ◀ Assembling the depth stop

The depth adjuster (CS 1324-2) with the locking nut (CS 1324-3) is screwed into the osmium depth stop (CS 1324-1). The drilling and tapping depth is determined by means of the x-ray images or CT/MRI images and adjusted on the depth stop. The depth stop is also used as a tissue protector. The drilling depth corresponds to the length of the screw.

CS 1217-1

CS 1210

CS 1217-3

CS 1324-1, -2, -3


CS 1323

CS 1325

# Surgical technique – Monocortical instrumentation

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## ◀ Bore for osmium screw

The drill hole is prepared using the drill bit, Ø 2.5 mm (CS 1323) through the osmium depth stop (CS 1324). Drill under C-arm control, if necessary.

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
## ◀ Tapping for osmium screw

The tap for the screw is prepared using the tap for osmium screws, Ø 5 mm (CS 1325) through the osmium depth stop.

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
## ◀ Preparing for screw insertion

First remove the bolt of the osmium screw with the screw driver (CS 1328). The attachment (CS 1310-2) is placed in the screw driver for osmium screw (CS 1310-1) and the osmium screw is locked. The prongs of the screw driver must be fitted to the grooves of the osmium screws.

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### ◀ Implanting the osmium screws

The osmium screw is implanted. The screw driver is removed and the second screw on the same side is inserted.

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### ◀ Implanting the osmium screws

After removal of the repositor and the reduction pins, the screws are inserted at the opposite side in the same way. Using the screw driver (CS 1310-1) or the tightener for osmium screws (CS 1322) the screws can be further tightened before being expanded, in order to achieve an even support of the wings.

CS 1217-1

CS 1210

CS 1310-1

CS 1310-2

CS 1310-1

CS 1322

# Surgical technique – Monocortical instrumentation

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## ◀ Fixing the osmium screws

The screw driver for osmium screw (CS 1310-1) without the attachment (CS 1310-2) is placed onto the osmium screw and the bolt is inserted into the shaft of the instrument. The bolt is locked using the screw driver, hex 2.5 mm (CS 1328) and the osmium screw is expanded.

**Important:** The screw driver for osmium screw must remain engaged on the grooves of the osmium screw, in order to prevent that the screw co-rotates.

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
## ◀ Fixation of the bone graft

If necessary the bone graft can be fixed using a cancellous bone screw in the whole in the center of the plate.

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# Surgical technique

## Surgical technique – Bicortical instrumentation

### Mounting the repositor, insertion of the plate, and compression

Use the same technique as with monocortical screws.

### Preparing the screw implantation

For the cancellous bone screws the same drill bit (CS 1322) and the depth stop (CS 1324) may be used. Determining the drill depth by using the depth stop is possible as well. The drilling depth corresponds to the screw length.

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### Tapping for cancellous screws

Use the same technique as with the monocortical screws. However, use the tap for cancellous bone screws, Ø 4.0 mm (CS 1326).

**Important:** It is optional to use the tap (CS 1326) before applying the screws. Use a C-arm.

### Implanting the cancellous bone screws


The appropriate length of the cancellous bone screws may be determined using the depth gauge (UT 1068-23). Inserting the cancellous bone screws by means of the screw driver (CS 1328) under C-arm control.

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







| Implants                                                   | Product number |
|------------------------------------------------------------|----------------|
| osmium screw incl. bolt, Ø 5 mm, length 14 mm, titanium    | CS 1300-14T    |
| osmium screw incl. bolt, Ø 5 mm, length 16 mm, titanium    | CS 1300-16T    |
| osmium screw incl. bolt, Ø 5 mm, length 18 mm, titanium    | CS 1300-18T    |
| Cancellous bone screw, Ø 4 mm, length 12 mm, titanium      | CS 1301-12T    |
| Cancellous bone screw, Ø 4 mm, length 14 mm, titanium      | CS 1301-14T    |
| Cancellous bone screw, Ø 4 mm, length 16 mm, titanium      | CS 1301-16T    |
| Cancellous bone screw, Ø 4 mm, length 18 mm, titanium      | CS 1301-18T    |
| Cancellous bone screw, Ø 4 mm, length 20 mm, titanium      | CS 1301-20T    |
| Cancellous bone screw, Ø 4 mm, length 22 mm, titanium      | CS 1301-22T    |
| Cancellous bone screw, Ø 4 mm, length 24 mm, titanium      | CS 1301-24T    |
| Cancellous bone screw, Ø 4 mm, length 26 mm, titanium      | CS 1301-26T    |
| osmium plate, 4 holes, length 23 mm, width 20 mm, titanium | CS 1306-23T    |
| osmium plate, 4 holes, length 27 mm, width 20 mm, titanium | CS 1306-27T    |
| osmium plate, 4 holes, length 30 mm, width 20 mm, titanium | CS 1306-30T    |
| osmium plate, 4 holes, length 33 mm, width 20 mm, titanium | CS 1306-33T    |
| osmium plate, 4 holes, length 36 mm, width 20 mm, titanium | CS 1306-36T    |
| osmium plate, 6 holes, length 40 mm, width 20 mm, titanium | CS 1307-40T    |
| osmium plate, 6 holes, length 44 mm, width 20 mm, titanium | CS 1307-44T    |
| osmium plate, 6 holes, length 48 mm, width 20 mm, titanium | CS 1307-48T    |
| osmium plate, 6 holes, length 52 mm, width 20 mm, titanium | CS 1307-52T    |
| osmium plate, 6 holes, length 56 mm, width 20 mm, titanium | CS 1307-56T    |
| osmium plate, 8 holes, length 60 mm, width 20 mm, titanium | CS 1308-60T    |
| osmium plate, 8 holes, length 66 mm, width 20 mm, titanium | CS 1308-66T    |
| osmium plate, 8 holes, length 72 mm, width 20 mm, titanium | CS 1308-72T    |

| Instruments                                           | Product number |
|-------------------------------------------------------|----------------|
| Reduction pin (two are necessary)                     | CS 1210        |
| Drill bit for reduction pin, Ø 2.0 mm                 | CS 1215        |
| Repositor                                             | CS 1217-1      |
| Drill guides for repositor (two are necessary)        | CS 1217-2      |
| Lever arm for repositor                               | CS 1217-3      |
| Screw holding forceps                                 | CS 1219        |
| Screw driver for osmium screw                         | CS 1310-1      |
| Attachment for screw driver                           | CS 1310-2      |
| Bone graft holder                                     | CS 1320        |
| Bone graft impactor                                   | CS 1321        |
| osmium screw tightener                                | CS 1322        |
| osmium drill bit, Ø 2.5 mm, for CS 1324               | CS 1323        |
| osmium depth stop                                     | CS 1324-1      |
| Depth adjustor for CS 1324-1                          | CS 1324-2      |
| Locking nut for CS 1324-1                             | CS 1324-3      |
| Tap for osmium screws, Ø 5 mm, for CS 1324            | CS 1325        |
| Tap for cancellous bone screws, Ø 4.0 mm, for CS 1324 | CS 1326        |
| Plate bender                                          | CS 1327-1      |
| Plate hold and bend instrument                        | CS 1327-2      |
| Screw driver, hex 2.5 mm                              | CS 1328        |
| Measuring caliper                                     | CS 5788        |
| Depth gauge, length 230 mm                            | UT 1068-23     |

**Ulrich GmbH & Co. KG**

Buchbrunnenweg 12

89081 Ulm

Germany

Telefon +49 (0)731 9654-110/225

Telefax +49 (0)731 9654-2702

e-mail [spine@ulrich-ulm.de](mailto:spine@ulrich-ulm.de)

internet [www.ulrich-ulm.de](http://www.ulrich-ulm.de)

WS 1301/10-03 Druckfehler und Irrtum vorbehalten · With reservation of misprints and errors · © Design Agentur Berstecher & Schoblocher · [www.design-agentur.de](http://www.design-agentur.de)