

Large Fragment

Large Fragment

4.01 - 4.14



You miss 100% of the
shots you don't take.



04

4.5mm Cortical Screws

Drill Bit : 3.2mm



4.5mm Cortical Screws, Self Tapping

Drill Bit : 3.2mm



04

Large Fragment

Length (mm)	Stainless Steel	Titanium
12	105.012	TI-105.012
14	105.014	TI-105.014
16	105.016	TI-105.016
18	105.018	TI-105.018
20	105.020	TI-105.020
22	105.022	TI-105.022
24	105.024	TI-105.024
26	105.026	TI-105.026
28	105.028	TI-105.028
30	105.030	TI-105.030
32	105.032	TI-105.032
34	105.034	TI-105.034
36	105.036	TI-105.036
38	105.038	TI-105.038
40	105.040	TI-105.040
42	105.042	TI-105.042
44	105.044	TI-105.044
46	105.046	TI-105.046
48	105.048	TI-105.048
50	105.050	TI-105.050
52	105.052	TI-105.052
54	105.054	TI-105.054
56	105.056	TI-105.056
58	105.058	TI-105.058
60	105.060	TI-105.060
62	105.062	TI-105.062
64	105.064	TI-105.064
66	105.066	TI-105.066
68	105.068	TI-105.068
70	105.070	TI-105.070
72	105.072	TI-105.072
74	105.074	TI-105.074
76	105.076	TI-105.076
78	105.078	TI-105.078
80	105.080	TI-105.080

Length (mm)	Stainless Steel	Titanium
12	105.212	TI-105.212
14	105.214	TI-105.214
16	105.216	TI-105.216
18	105.218	TI-105.218
20	105.220	TI-105.220
22	105.222	TI-105.222
24	105.224	TI-105.224
26	105.226	TI-105.226
28	105.228	TI-105.228
30	105.230	TI-105.230
32	105.232	TI-105.232
34	105.234	TI-105.234
36	105.236	TI-105.236
38	105.238	TI-105.238
40	105.240	TI-105.240
42	105.242	TI-105.242
44	105.244	TI-105.244
46	105.246	TI-105.246
48	105.248	TI-105.248
50	105.250	TI-105.250
52	105.252	TI-105.252
54	105.254	TI-105.254
56	105.256	TI-105.256
58	105.258	TI-105.258
60	105.260	TI-105.260
62	105.262	TI-105.262
64	105.264	TI-105.264
66	105.266	TI-105.266
68	105.268	TI-105.268
70	105.270	TI-105.270
72	105.272	TI-105.272
74	105.274	TI-105.274
76	105.276	TI-105.276
78	105.278	TI-105.278
80	105.280	TI-105.280

4.5mm Malleolar Screws

Drill Bit : 3.2mm



Length (mm)	Stainless Steel	Titanium
25	106.025	TI-106.025
30	106.030	TI-106.030
35	106.035	TI-106.035
40	106.040	TI-106.040
45	106.045	TI-106.045
50	106.050	TI-106.050
55	106.055	TI-106.055
60	106.060	TI-106.060
65	105.065	TI-105.065
70	106.070	TI-106.070
75	106.075	TI-106.075
80	106.080	TI-106.080
85	106.085	TI-106.085
90	106.090	TI-106.090
95	106.095	TI-106.095
100	106.100	TI-106.100

Washers for 4.5mm to 7mm Screws



	Stainless Steel	Titanium
Code	100.120	TI-100.120

6.5mm Cancellous Screws, 16mm Thread

Drill Bit : 3.2mm



Length (mm)	Stainless Steel	Titanium
25	111.025	TI-111.025
30	111.030	TI-111.030
35	111.035	TI-111.035
40	111.040	TI-111.040
45	111.045	TI-111.045
50	111.050	TI-111.050
55	111.055	TI-111.055
60	111.060	TI-111.060
65	111.065	TI-111.065
70	111.070	TI-111.070
75	111.075	TI-111.075
80	111.080	TI-111.080
85	111.085	TI-111.085
90	111.090	TI-111.090
95	111.095	TI-111.095
100	111.100	TI-111.100
105	111.105	TI-111.105
110	111.110	TI-111.110
115	111.115	TI-111.115
120	111.120	TI-111.120

6.5mm Cancellous Screws, 32mm Thread

Drill Bit : 3.2mm



Length (mm)	Stainless Steel	Titanium
40	112.040	TI-112.040
45	112.045	TI-112.045
50	112.050	TI-112.050
55	112.055	TI-112.055
60	112.060	TI-112.060
65	112.065	TI-112.065
70	112.070	TI-112.070
75	112.075	TI-112.075
80	112.080	TI-112.080
85	112.085	TI-112.085
90	112.090	TI-112.090
95	112.095	TI-112.095
100	112.100	TI-112.100
105	112.105	TI-112.105
110	112.110	TI-112.110
115	112.115	TI-112.115
120	112.120	TI-112.120

Large Fragment - Bone Screws



6.5mm Cancellous Screws, Full Thread

Drill Bit : 3.2mm



Length (mm)	Stainless Steel	Titanium
25	113.025	TI-113.025
30	113.030	TI-113.030
35	113.035	TI-113.035
40	113.040	TI-113.040
45	113.045	TI-113.045
50	113.050	TI-113.050
55	113.055	TI-113.055
60	113.060	TI-113.060
65	113.065	TI-113.065
70	113.070	TI-113.070
75	113.075	TI-113.075
80	113.080	TI-113.080
85	113.085	TI-113.085
90	113.090	TI-113.090
95	113.095	TI-113.095
100	113.100	TI-113.100
105	113.105	TI-113.105
110	113.110	TI-113.110
115	113.115	TI-113.115
120	113.120	TI-113.120

IS003 Screw Set Large

Item Codes	Set Consisting of	Units
105.014-105.068	4.5mm Cortical Screw 6 each (Length in mm) - 14, 16, 18, 20, 22, 24 7 each (Length in mm) - 26, 28, 38, 42, 44, 46 14 each (Length in mm)- 30, 32, 34, 36, 40 3 each (Length in mm) - 48, 52, 56, 60, 64, 68	166
106.025-106.070	4.5mm Mallelor Screws 3 each (Length in mm) - 25, 30, 35, 40, 45, 50, 55, 60, 65, 70	30
111.030-111.110	6.5mm Cancellous Bone Screw, 16mm Thread 2 each (Length in mm) - 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110	34
112.045-112.110	6.5mm Cancellous Bone Screw, 32mm Thread 3 each (Length in mm)- 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110	42
113.025-113.110	6.5mm Cancellous Bone Screw, Fully Threaded 4 each (Length in mm)- 25, 30, 35, 40, 2 each (Length in mm)- 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110	44
100.120	Washer for Large Screw	12
3415.000	Screw Holding Forceps	1
A10.006	Aluminium Box for Large Screw Set	1

4.5mm Narrow Dynamic Compression Plate

Fixation with 4.5 & 6.5mm screws.
All plates are designed as self compression holes.
The end holes are designed for cancellous bone screws.



Holes	Stainless Steel	Titanium
4	722.104	TI-722.104
5	722.105	TI-722.105
6	722.106	TI-722.106
7	722.107	TI-722.107
8	722.108	TI-722.108
9	722.109	TI-722.109
10	722.110	TI-722.110
11	722.111	TI-722.111
12	722.112	TI-722.112
13	722.113	TI-722.113
14	722.114	TI-722.114
15	722.115	TI-722.115
16	722.116	TI-722.116
17	722.117	TI-722.117
18	722.118	TI-722.118

4.5mm Narrow Limited Contact Dynamic Compression Plate

Fixation with 4.5 & 6.5mm screws.
The LC-DCP holes allow compression in both the direction of the longitudinal axis.
The end holes are suitable for cancellous bone screws.
Undercuts for vascular preservation



Holes	Stainless Steel	Titanium
4	723.104	TI-723.104
5	723.105	TI-723.105
6	723.106	TI-723.106
7	723.107	TI-723.107
8	723.108	TI-723.108
9	723.109	TI-723.109
10	723.110	TI-723.110
11	723.111	TI-723.111
12	723.112	TI-723.112
13	723.113	TI-723.113
14	723.114	TI-723.114
15	723.115	TI-723.115
16	723.116	TI-723.116
17	723.117	TI-723.117
18	723.118	TI-723.118

4.5mm Broad Dynamic Compression Plate


Fixation with 4.5 & 6.5mm Screws.
All Plates are designed as self compression holes.
The end holes are suitable for cancellous bone screws.



Holes	Stainless Steel	Titanium
4	724.104	TI-724.104
5	724.105	TI-724.105
6	724.106	TI-724.106
7	724.107	TI-724.107
8	724.108	TI-724.108
9	724.109	TI-724.109
10	724.110	TI-724.110
11	724.111	TI-724.111
12	724.112	TI-724.112
13	724.113	TI-724.113
14	724.114	TI-724.114
15	724.115	TI-724.115
16	724.116	TI-724.116
17	724.117	TI-724.117
18	724.118	TI-724.118
19	724.119	TI-724.119
20	724.120	TI-724.120

4.5mm Broad Limited Contact Dynamic Compression Plate

Fixation with 4.5 & 6.5mm Screws.
The LC-DCP holes allow compression in both the direction of the longitudinal axis.
The end holes are suitable for cancellous bone screws.
Undercuts for vascular preservation.



Holes	Stainless Steel	Titanium
4	725.104	TI-725.104
5	725.105	TI-725.105
6	725.106	TI-725.106
7	725.107	TI-725.107
8	725.108	TI-725.108
9	725.109	TI-725.109
10	725.110	TI-725.110
11	725.111	TI-725.111
12	725.112	TI-725.112
13	725.113	TI-725.113
14	725.114	TI-725.114
15	725.115	TI-725.115
16	725.116	TI-725.116
17	725.117	TI-725.117
18	725.118	TI-725.118
19	725.119	TI-725.119
20	725.120	TI-725.120

4.5mm Narrow Lengthening Plate, 8 Holes



Lengthening (mm)	Length (mm)	Stainless Steel	Titanium
30	135	755.030	TI-755.030
40	145	755.040	TI-755.040
50	155	755.050	TI-755.050
60	165	755.060	TI-755.060
70	175	755.070	TI-755.070
80	185	755.080	TI-755.080

4.5mm Broad Lengthening Plate, 8 Holes



Lengthening (mm)	Length (mm)	Stainless Steel	Titanium
30	135	754.030	TI-754.030
40	145	754.040	TI-754.040
50	155	754.050	TI-754.050
60	165	754.060	TI-754.060

4.5mm Broad Lengthening Plate, 10 Holes



Lengthening (mm)	Length (mm)	Stainless Steel	Titanium
50	179	753.050	TI-753.050
60	189	753.060	TI-753.060
70	199	753.070	TI-753.070
80	209	753.080	TI-753.080
90	219	753.090	TI-753.090
100	229	753.100	TI-753.100
110	239	753.110	TI-753.110
120	249	753.120	TI-753.120

4.5mm Semi Tubular Plates



Holes	Stainless Steel	Titanium
4	704.204	TI-704.204
5	704.205	TI-704.205
6	704.206	TI-704.206
7	704.207	TI-704.207
8	704.208	TI-704.208
9	704.209	TI-704.209
10	704.210	TI-704.210
11	704.211	TI-704.211
12	704.212	TI-704.212

4.5mm Lateral Tibial Head Buttress Plate for Left Leg

Fixation with 4.5 & 6.5mm Screws



Holes	Stainless Steel	Titanium
5	733.105R	TI-733.105R
7	733.107R	TI-733.107R
9	733.109R	TI-733.109R
11	733.111R	TI-733.111R
13	733.113R	TI-733.113R

4.5mm Lateral Tibial Head Buttress Plate for Right Leg

Fixation with 4.5 & 6.5mm Screws



Holes	Stainless Steel	Titanium
5	733.105L	TI-733.105L
7	733.107L	TI-733.107L
9	733.109L	TI-733.109L
11	733.111L	TI-733.111L
13	733.113L	TI-733.113L

4.5mm Condylar Buttress Plate, Left



Holes	Stainless Steel	Titanium
7	731.107L	TI-731.107L
9	731.109L	TI-731.109L
11	731.111L	TI-731.111L
12	731.112L	TI-731.112L
13	731.113L	TI-731.113L
15	731.115L	TI-731.115L

4.5mm Condylar Buttress Plate, Right



Holes	Stainless Steel	Titanium
7	731.107R	TI-731.107R
9	731.109R	TI-731.109R
11	731.111R	TI-731.111R
12	731.112R	TI-731.112R
13	731.113R	TI-731.113R
15	731.115R	TI-731.115R

4.5mm T-Plate



Holes	Stainless Steel	Titanium
3	713.103	TI-713.103
4	713.104	TI-713.104
5	713.105	TI-713.105
6	713.106	TI-713.106
7	713.107	TI-713.107
8	713.108	TI-713.108
9	713.109	TI-713.109
10	713.110	TI-713.110
11	713.111	TI-713.111
12	713.112	TI-713.112
13	713.113	TI-713.113
14	713.114	TI-713.114

4.5mm T-Plate Buttress



Holes	Stainless Steel	Titanium
4	714.104	TI-714.104
5	714.105	TI-714.105
6	714.106	TI-714.106
7	714.107	TI-714.107
8	714.108	TI-714.108
9	714.109	TI-714.109
10	714.110	TI-714.110
11	714.111	TI-714.111
12	714.112	TI-714.112
13	714.113	TI-714.113
14	714.114	TI-714.114

4.5mm L-Plate Buttress, Right Angled
(for Left Leg)



Holes	Stainless Steel	Titanium
3	706.103R	TI-706.103R
4	706.104R	TI-706.104R
5	706.105R	TI-706.105R
6	706.106R	TI-706.106R
7	706.107R	TI-706.107R
8	706.108R	TI-706.108R
9	706.109R	TI-706.109R
10	706.110R	TI-706.110R
11	706.111R	TI-706.111R
12	706.112R	TI-706.112R
13	706.113R	TI-706.113R
14	706.114R	TI-706.114R

4.5 mm L-Plate Buttress, Left Angled
(for Right Leg)



Holes	Stainless Steel	Titanium
3	706.103L	TI-706.103L
4	706.104L	TI-706.104L
5	706.105L	TI-706.105L
6	706.106L	TI-706.106L
7	706.107L	TI-706.107L
8	706.108L	TI-706.108L
9	706.109L	TI-706.109L
10	706.110L	TI-706.110L
11	706.111L	TI-706.111L
12	706.112L	TI-706.112L
13	706.113L	TI-706.113L
14	706.114L	TI-706.114L

4.5mm Spoon Plate

Fixation with 4.5mm Cortical Screws
Shaft with V-profile for the fixation on the frontal edge of the tibia.



Holes	Stainless Steel	Titanium
5	718.105	TI-718.105
6	718.106	TI-718.106

4.5mm Reconstruction Plate, Straight

Do not bend more than 15°



Holes	Stainless Steel	Titanium
3	728.103	TI-728.103
4	728.104	TI-728.104
5	728.105	TI-728.105
6	728.106	TI-728.106
7	728.107	TI-728.107
8	728.108	TI-728.108
9	728.109	TI-728.109
10	728.110	TI-728.110
11	728.111	TI-728.111
12	728.112	TI-728.112
13	728.113	TI-728.113
14	728.114	TI-728.114
15	728.115	TI-728.115
16	728.116	TI-728.116
17	728.117	TI-728.117
18	728.118	TI-728.118
20	728.120	TI-728.120
22	728.122	TI-728.122

4.5mm Cobra Head Plate



Holes	Stainless Steel	Titanium
8	735.108	TI-735.108
9	735.109	TI-735.109
10	735.110	TI-735.110
11	735.111	TI-735.111
12	735.112	TI-735.112
13	735.113	TI-735.113
14	735.114	TI-735.114

4.5mm Hook Plate

Fixation with 4.5 & 6.5mm Screws
Shaft with DCP dynamic compression holes and slot for tension device. Used after a blade of an angled hip plate breaks off at the proximal femur.



Holes	Stainless Steel	Titanium
6	730.106	TI-730.106

Roy Camille Plate



Holes	Stainless Steel	Titanium
5	804.005	TI-804.005
6	804.006	TI-804.006
7	804.007	TI-804.007
8	804.008	TI-804.008
9	804.009	TI-804.009
10	804.010	TI-804.010
11	804.011	TI-804.011
12	804.012	TI-804.012
13	804.013	TI-804.013
14	804.014	TI-804.014

Sherman Plate



Holes	Stainless Steel	Titanium
3	802.003	TI-802.003
4	802.004	TI-802.004
5	802.005	TI-802.005
6	802.006	TI-802.006
7	802.007	TI-802.007
8	802.008	TI-802.008
9	802.009	TI-802.009
10	802.010	TI-802.010
11	802.011	TI-802.011
12	802.012	TI-802.012
13	802.013	TI-802.013
14	802.014	TI-802.014

4.5mm Proximal Plate Lateral Side



Holes	Type	Stainless Steel	Titanium
3	Right	MP04.03R	TI-MP04.03R
4	Right	MP04.04R	TI-MP04.04R
5	Right	MP04.05R	TI-MP04.05R
3	Left	MP04.03L	TI-MP04.03L
4	Left	MP04.04L	TI-MP04.04L
5	Left	MP04.05L	TI-MP04.05L

4.5mm Proximal Tibia Plate with Round Holes



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
4	MP03.04R	TI-MP03.04R	MP03.04L	TI-MP03.04L
5	MP03.05R	TI-MP03.05R	MP03.05L	TI-MP03.05L
6	MP03.06R	TI-MP03.06R	MP03.06L	TI-MP03.06L
7	MP03.07R	TI-MP03.07R	MP03.07L	TI-MP03.07L
8	MP03.08R	TI-MP03.08R	MP03.08L	TI-MP03.08L
9	MP03.09R	TI-MP03.09R	MP03.09L	TI-MP03.09L
10	MP03.10R	TI-MP03.10R	MP03.10L	TI-MP03.10L

4.5mm Proximal Plate Medial Side



Holes	Type	Stainless Steel	Titanium
4	Right	MP01.04R	TI-MP01.04R
5	Right	MP01.05R	TI-MP01.05R
6	Right	MP01.06R	TI-MP01.06R
7	Right	MP01.07R	TI-MP01.07R
4	Left	MP01.04L	TI-MP01.04L
5	Left	MP01.05L	TI-MP01.05L
6	Left	MP01.06L	TI-MP01.06L
7	Left	MP01.07L	TI-MP01.07L

4.5mm Proximal Lateral Tibial Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	736.103L	TI-736.103L	736.103R	TI-736.103R
5	736.105L	TI-736.105L	736.105R	TI-736.105R
7	736.107L	TI-736.107L	736.107R	TI-736.107R
9	736.109L	TI-736.109L	736.109R	TI-736.109R
11	736.111L	TI-736.111L	736.111R	TI-736.111R
13	736.113L	TI-736.113L	736.113R	TI-736.113R

4.5mm Proximal Tibia Medial Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	763.103L	TI-763.103L	763.103R	TI-763.103R
5	763.105L	TI-763.105L	763.105R	TI-763.105R
7	763.107L	TI-763.107L	763.107R	TI-763.107R
9	763.109L	TI-763.109L	763.109R	TI-763.109R
11	763.111L	TI-763.111L	763.111R	TI-763.111R
13	763.113L	TI-763.113L	763.113R	TI-763.113R

4.5mm Proximal Femur Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	764.103L	TI-764.103L	764.103R	TI-764.103R
5	764.105L	TI-764.105L	764.105R	TI-764.105R
7	764.107L	TI-764.107L	764.107R	TI-764.107R
9	764.109L	TI-764.109L	764.109R	TI-764.109R
11	764.111L	TI-764.111L	764.111R	TI-764.111R
13	764.113L	TI-764.113L	764.113R	TI-764.113R

4.5mm Distal Tibial Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	737.103L	TI-737.103L	737.103R	TI-737.103R
5	737.105L	TI-737.105L	737.105R	TI-737.105R
7	737.107L	TI-737.107L	737.107R	TI-737.107R
9	737.109L	TI-737.109L	737.109R	TI-737.109R
11	737.111L	TI-737.111L	737.111R	TI-737.111R
13	737.113L	TI-737.113L	737.113R	TI-737.113R

4.5mm Distal Tibial Lateral Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	737.203L	TI-737.203L	737.203R	TI-737.203R
5	737.205L	TI-737.205L	737.205R	TI-737.205R
7	737.207L	TI-737.207L	737.207R	TI-737.207R
9	737.209L	TI-737.209L	737.209R	TI-737.209R
11	737.211L	TI-737.211L	737.211R	TI-737.211R
13	737.213L	TI-737.213L	737.213R	TI-737.213R

4.5mm Distal Lateral Tibia Plate with Round Holes



Holes	Type	Stainless Steel	Titanium
5	Left	CP01.05L	TI-CP01.05L
5	Right	CP01.05R	TI-CP01.05R
7	Left	CP01.07L	TI-CP01.07L
7	Right	CP01.07R	TI-CP01.07R
9	Left	CP01.09L	TI-CP01.09L
9	Right	CP01.09R	TI-CP01.09R
11	Left	CP01.11L	TI-CP01.11L
11	Right	CP01.11R	TI-CP01.11R
13	Left	CP01.13L	TI-CP01.13L
13	Right	CP01.13R	TI-CP01.13R

4.5mm Distal Femur Plate



Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	737.303L	TI-737.303L	737.303R	TI-737.303R
5	737.305L	TI-737.305L	737.305R	TI-737.305R
7	737.307L	TI-737.307L	737.307R	TI-737.307R
9	737.309L	TI-737.309L	737.309R	TI-737.309R
11	737.311L	TI-737.311L	737.311R	TI-737.311R
13	737.313L	TI-737.313L	737.313R	TI-737.313R

4.5mm Distal Humeral Plate



Holes	Type	Stainless Steel	Titanium
4	Right	MP02.04R	TI-MP02.04R
5	Right	MP02.05R	TI-MP02.05R
6	Right	MP02.06R	TI-MP02.06R
9	Right	MP02.09R	TI-MP02.09R
4	Left	MP02.04L	TI-MP02.04L
5	Left	MP02.05L	TI-MP02.05L
6	Left	MP02.06L	TI-MP02.06L
9	Left	MP02.09L	TI-MP02.09L

4.5mm Fibular Distal Tibia Plate



Holes	Type	Stainless Steel	Titanium
7	Right	CP12.07R	TI-CP12.07R
9	Right	CP12.09R	TI-CP12.09R
11	Right	CP12.11R	TI-CP12.11R
13	Right	CP12.13R	TI-CP12.13R
7	Left	CP12.07L	TI-CP12.07L
9	Left	CP12.09L	TI-CP12.09L
11	Left	CP12.11L	TI-CP12.11L
13	Left	CP12.13L	TI-CP12.13L